BT202

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BT202

Self-contained Biological Indicator.

Rapid Readout Colorimetric System.

Usage

Monitoring Vacuum-assisted and Gravity air displacement Steam Sterilization Processes at 121-135 °C.

Applicable Regulation

ISO 11138-1:2006 and ISO 11138-3:2006; IRAM 37102-1:1999 and IRAM 37102-3:1999.

Classification

Class 1, according to risk.

Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

Characteristics

Polypropylene random copolymer tube: 48.0 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm.

Polypropylene cap: 10.5 mm diameter, 13.0 mm high. Wall thickness: 1.0 mm.

Glass ampoule: 42.0 mm high. mm. External diameter: 6.5 mm. Wall thickness: 0.2 – 0.3 mm.

Culture medium: 0.7 ml, **purple** color.

Cap filter: medical grade paper, 60 g/m², 16.0 mm diameter.

10⁵ or 10⁶ *Geobacillus stearothermophilus* ATCC 7953 spores per vial.

Spores Carrier: Filter Paper, 8.0 mm diameter, 160 g/m².

Positive growth results (failure in sterilization process): Culture medium turns to **yellow**, incubation at 60 °C.

Sterilization failures might be detected within 3 to 5 hours of incubation (sensitivity ≥ 97 %). Recommended Incubation Time: 10 hours.

D-Value: Not lower than 1.5 minutes at 121 °C.

Another D-value is estimated at 134 °C for the calculation of Z-Value

Z-Value: Not lower than 6 °C.

Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place.

Transportation conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

Shelf-life

24 months

Packing

100 units per box.

Packing information: product code and description, process for intended use, presentation, classification according to regulation, manufacturer information and data on box's label. Note: manufacture date is calculated by subtracting 24 months to the expiration date.

Weight per box: 241,2 g.

Labelling

On product: 15.0 mm x 33.0 mm polypropylene label. Printed in black. 1.5 mm chemical indicator line, printed with Steam reactive ink (Color changes to brown). Product code and batch number, process for intended use and bacterial strain.

On product's box: product code, batch number, bacterial load, manufacture and expiration date.

Possible target markets

Healthcare and Industry.

Other important information (optional)

Read product's directions for use thoroughly before use.

Note: When necessary, limits of each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.